

PRICE-LIST

w.e.f JAN 2022





PRECISION BIOMED PVT. LTD. an ISO 9001:2015/ ISO 13485:2016/CE & WHO GMP Certified company with a flair for Creativity & Professionalism.

We are well equipped with modern state of Art & Technology dealing in Diagnostic products and Equipments. Our Manufacturing Plant is situated at Jaipur which is nearby Delhi with all modern facility in Approx 20000 sq. ft. area as per new medical device Rules 2017. We are doing manufacturing of all Biochemistry/Serology and Rapid Test Kits.

PRECISION BIOMED PVT. LTD. has started its operation in the year 2002 with **Micropipettes** and are expanding with **Serology, Urine Strips, Biochemistry, Rapid Tests and Instruments.**

PRECISION BIOMED PVT. LTD. is a young company where its people combines the best of youth and maturity. Due to its diversified body structure it has earned a respectable place in Diagnostic industry We are proud to say that **Precision** is the only company to give **HIGH CLASS QUALITY PRODUCTS** on Economical price.

All departments' i.e. Q.C. Servicing and Production work in perfect harmony to fulfill expectation and standards of our worthy customers.

Value of money and sharing our customer's vision remain the core of our working Philosophy.

BIOCHEMISTRY (LIQUID)

S.No.	Products	Pack Size	Kit MRP
1.	Albumin	2x50 ML	482.00
2.	Albumin	2x500 ML	3,614.00
3.	Amylase	2x10 ML	2,477.00
4.	Amylase	2x25 ML	4,954.00
5.	Amylase	2x50 ML	9,908.00
6.	Alkaline Phosphatase	2x25 ML	833.00
7.	Alkaline Phosphatase	2x50 ML	1,515.00
8.	Bilirubin (Total & Direct)	2 x 100 ML	895.00
9.	Bilirubin (Total & Direct)	2x500 ML	3,579.00
10.	Cholesterol	2x25 ML	719.00
11.	Cholesterol	2x50 ML	1,308.00
12.	Cholesterol	5x100 ML	5,230.00
13.	Cholesterol	2x500 ML	9,414.00
14.	CK-NAC	10 ML	1,287.00
15.	CK-NAC	25 ML	2,926.00
16.	CK-MB	10 ML	1,969.00
17.	CK-MB	25 ML	4,474.00
18.	Creatinine	2X50 ML	895.00
19.	Creatinine	2x500 ML	6,710.00
20.	Creatinine-Single	2x50 ML	895.00
21.	Creatinine-Single	2x500 ML	6,710.00





BIOCHEMISTRY (LIQUID)

S.No.	Products	Pack Size	Kit MRP
22.	Gamma-GT	10 ML	633.00
23.	Gamma-GT	25 ML	1,438.00
24.	Gamma-GT	5x100ML	14,380.00
25.	Glucose	5x100 ML	946.00
26.	Glucose	2x500 ML	1,703.00
27.	LDH	10 ML	759.00
28.	LDH	25 ML	1,725.00
29.	LDH	5x100 ML	17,250.00
30.	SGOT (AST)	2x25 ML	833.00
31.	SGOT (AST)	2x50 ML	1,515.00
32.	SGOT (AST)	5x100 ML	6058.00
33.	SGOT (AST)	2x500 ML	10,904.00
34.	SGPT (ALT)	2x25 ML	833.00
35.	SGPT (ALT)	2x50ML	1,515.00
36.	SGPT (ALT)	5x100ML	6,058.00
37.	SGPT (ALT)	2x500 ML	10,904.00
38.	Triglycerides	2x25 ML	1,514.00
39.	Triglycerides	2x50 ML	2,753.00
40.	Triglycerides	5x100 ML	11,012.00
41.	Triglycerides	2x500ML	19,822.00
42.	Total Protein	2x50 ML	482.00
43.	Total Protein	2x500 ML	3,614.00
44.	Urea Berthelot	2x50 ML	895.00
45.	Urea Berthelot	2x500 ML	6,710.00
46.	Urea-UV	2x25 ML	833.00
47.	Urea-UV	2x50 ML	1,515.00
48.	Urea-UV	5x100 ML	6,058.00
49.	Urea-UV	2x500 ML	10,904.00
50.	URIC ACID	2x25 ML	1,136.00
51.	URIC ACID	2x50 ML	2,065.00
52.	URIC ACID	5x100 ML	9,292.00
53.	URIC ACID	2x500 ML	16,726.00

TURBI LATEX

1.	TURBI-ASO	1x50 ML	6,584.00
2.	TURBI-CRP	1X50 ML	5,985.00
3.	TURBI-RF	1X50 ML	5,387.00



MONO-VIAL

S.No.	Products	Pack Size	Kit MRP
1.	SODIUM - MONO	50X1 ML	1,584.00
2.	POTASSIUM-MONO	50X1 ML	1,584.00
3.	PHOSPHORUS-MONO	50 X 1 ML	1,584.00
4.	MAGNESIUM-MONO	50X1 ML	2,477.00
5.	CALCIUM-MONO	50X1 ML	1,033.00
6.	LACTATE-MONO	25X1 ML	7,728.00
7.	CHLORIDE-MONO	50X1 ML	1,480.00
8.	ALK. PHOSPHATE-MONO	50X1 ML	1,584.00
9.	AMYLASE-MONO	25X1 ML	3,400.00

SPECIALIZED PRODUCTS

1.	HDL CHOLESTEROL - DIRECT	20 ML	2,753.00
2.	HDL CHOLESTEROL - DIRECT	40 ML	4,955.00
3.	HDL CHOLESTEROL-DIRECT	80 ML	8,810.00
4.	HDL CHOLESTEROL-DIRECT (SINGLE)	20 ML	2,850.00
5.	HDL CHOLESTEROL-DIRECT (SINGLE)	40 ML	5,130.00
6.	HDL CHOLESTEROL-DIRECT (SINGLE)	80 ML	9,066.00
7.	LDL CHOLESTEROL-DIRECT	20 ML	3,441.00
8.	LDL CHOLESTEROL-DIRECT	40 ML	6,194.00
9.	LDL CHOLESTEROL-DIRECT	80 ML	11,011.00
10.	LIPASE	10 ML	2,974.00
11.	LIPASE	20 ML	5,948.00
12.	ADA-MTB (KINETIC)	15 ML	5,888.00
13.	ADA-MTB (KINETIC)	30 ML	10,706.00
14.	ADA WITH CALIBRATOR	30 ML	11,868.00
15.	ADA WITH CALIBRATOR	90 ML	32,120.00
16.	AMMONIA-ULTRA	20 ML	5,506.00

COUGULATION

1.	PROTHROMBIN TIME	1X5 ML	778.00
2.	PROTHROMBIN TIME	4X5 ML	2,633.00

SEROLOGY

1.	ASO	25 Test	1,077.00
2.	ASO	100 Test	3,591.00
3.	CRP	25 Test	1,017.00



SEROLOGY

S.No.	Products	Pack Size	Kit MRP
4.	CRP	100 Test	2,633.00
5.	RF	25 Test	539.00
6.	RF	100 Test	1,915.00
7.	WIDAL SLIDE	4 x 5 ml	479.00
8.	WIDAL SLIDE ("O" & "H")	2X2X5 ML	539.00

BLOOD GROUPING

1.	ANTI A	10 ML	200.00
2.	ANTI B	10 ML	200.00
3.	ANTI D (IgM)	10 ML	450.00
4.	ANTI ABD (IgM)	3x10 ML	850.00

RAPID CARD TEST

1.	MALARIA PF/PV ANTIGEN CARD	30 Test	2,400.00
2.	DENGUE IgG/IgM	10 Test	2,519.00
3.	DENGUE IgG/IgM	25 Test	5,668.00
4.	DENGUE Ns1	10 Test	3,500.00
5.	DENGUE Ns1	25 Test	7,875.00
6.	DENGUE COMBO (IgG/IgM/NS1)	10 Test	5,500.00
7.	DENGUE COMBO (IgG/IgM/NS1)	25 Test	12,375.00
8.	HBSAG	50 Test	1,750.00
9.	HCV CARD	50 Test	4,500.00
10.	HIV 1+2 TRI-LINE CARD	50 Test	4,500.00
11.	TYPHOID IgM/IgM CARD	10 Test	800.00
12.	TYPHOID IgM/IgM CARD	25 Test	1,800.00
13.	URINE STRIP -GP	100 STRIPS	425.00
14.	CHIKUNGUNYA CARD	10 Test	5,700.00
15.	CHIKUNGUNYA CARD	25 Test	12,825.00
16.	SYPHILIS CARD	30 Test	900.00





1.	ECO Micropipettes (Fixed)	Each	-
2.	ECO Micropipettes (Variable)	Each	-
3.	PRO Micropipettes (Fixed)	Each	-
4.	PRO Micropipettes (Variable)	Each	-
5.	CLASSI Micropipettes (Fixed)	Each	-
6.	CLASSI Micropipettes (Variable)	Each	-
7.	AUTOCLAVABLE Micropipettes (Fixed)	Each	-
8.	AUTOCLAVABLE Micropipettes (Variable)	Each	-
9.	Mini Pipette	Each	-

INSTRUMENTS

S.No.	Products	Pack Size	Kit MRP
1.	PRECICHEM - 300+ (Semi Auto Analyzer)	Unit	-
2.	PRECINEPHCHEM (Protein Analyzer)	Unit	-
3.	PRECINEPHCHEM + (Protein Analyzer)	Unit	-
4.	PRECIELISA (Elisa Reader)	Unit	-
5.	PRECILYTE-i + (Electrolyte Analyzer)	Unit	-
6.	PRECICELL-300 (Hematology Analyzer)	Unit	-
7.	Dry Bath Incubator	Unit	-



ANALYZERS



PreciCELL-300


- Fully Automatic Double Chamber
- 3-Part Differential Cell Counting
- 20 Parameters + 3 Histograms
- Throughput : 60 Samples Per Hour
- Sample Volume : 20 µl
- Manual & Auto Calibration
- Patient Results Up to 1000 Samples
- Built In Thermal Printer
- Cyanide free, Eco Friendly Reagent
- Easy to Use, User Friendly Operating Menu

SPECIFICATIONS

Analysis principle	: Ion selective electrode (ISE)
Analysis Time	: less than 60 seconds
Calibration	: Automatic or Manual, 2-point
Calibrators	: Calibrator 1 (480 ml) & Calibrator 2 (100 ml) Individual bottle
system	
Throughput	: 60 Tests/Hr
Temperature	: 100 C - 400 C, Humidity > 85%
Sample Type	: Whole Blood, Serum, Plasma
Sample Input	: From tube, syringe, capillary
Sample size	: 100µL
Measurement	: Direct "Flow Through" (ISE)
Rinse method	: Automatic
Patient details	: Name Patient ID (keyboard/barcode) or auto-increment
Printer	: Thermal 58mm built-in printer
Data storage	: Up to 20000
QC data capacity	: 1000 sample per QC
Data transfer	: On-screen recall, Pen drive (PDF), Direct to PC (word/txt file)
Power requirement	: 165-242V, 50/60Hz, Grounding is compulsory
Dimension	: 285 x 202 x 350mm
Weight	: 6 Kg



PreciLYTE-i+
(Electrolyte Analyzer)



Preci Nephchem
(A Specific Protein Analyzer/
Nephelo-Chemistry Analyzer)

Features & Benefits

- Easy Operation, and no need of multi-point calibration.
- Simple Two step procedure all specific Protein Parameters eases handle of Operation.
- Value of money, each test packs can be assured to finish as per the number of samples without any wastage.
- Quality Control validation can be available as additional request by customer.
- The Liner Range is wide with high Stability.
- Professionally designed function of hardware and software can eliminate interference with thermal fluctuations.
- Equipped with automatic mixing system, which ensure active mixing.
- The external incubators is well mentioned by desired optimum temperature of 37°C for both cuvettes and Reagents as Pre-warm of all Antigen-Antibody Chemistries

PreciCHEM-300+

(Semi Automated Clinical Chemistry Analyzer)

- Robust and reliable diagnostics solution
- Consistent performance and easy operation
- Internal Incubator and coagulation test (PT/APIT/TT)



SOFTWARE FEATURES

PreciCHEM software is designed to simplify daily laboratory routine. It's simple and intuitive work flow enhances user efficiency and productivity

INTUITIVE USER INTERFACE

- Six informative screens control operations
- Easy navigation through submenus
- Requires minimal training

QUALITY CONTROL

- Mean, SD, % CV and R calculated for each assay
- Levey-Jennings control chart
- 3 levels per assay
- Westgard Rules

CALIBRATION

- Factor, one-point, two-point and multi-point
- Non liner multi point calibration
- Up to 3 replicates per point
- Up to 10 calibration points

ERROR DETECTION

- Early detection and reporting of errors on screen
- Saves time and avoids faulty test runs which saves reagents

LIGHT SOURCE
Halogen Quartz tungsten lamp (6V, 10Watt)
Lamp life saving mode (programmable sleep interval)
Unique Lamp life enhancing circuit

ASPIRATION
High precision peristaltic pump
Non corrosive tubing ensure long life and accurate aspiration
Aspiration air gap and auto empty function after each sample run
Minimal chances of carry over

EASY MAINTENANCE
Self-diagnosis and system health reporting are an integral part of PreciLab. Easily set photometer gain and perform pump calibration

ONE TOUCH ACCESS TO DEDICATED ASSAYS
48 direct access assays

PRINTABLE REAL TIME GRAPH
Real time graph is display on screen and can be printed

TERMS & CONDITIONS

Price	:	The above prices are effective from 1st January 2022 and are subject to change without notice.
Supply	:	Supplies will be made at prices prevailing at the time of dispatch, All supplies will be made on terms cheque against delivery or through bank.
Taxes	:	The prices mentioned in these price list are exclusive of GST or any other levies, which will be charged extra, as applicable
Products	:	All products are checked and packed with utmost care. No claim or breakage, delay or damage will be entertained after the items have left our premises. Goods once sold will not be accepted for credit or exchange .
Storage	:	Products have to be stored at temperatures Indicated on each pack.
Delivery	:	Immediate if in stock or 4 to 5 weeks after receipt of confirmed order.
Transport	:	Packaging, Forwarding, Postage, Transportation & Incidental charges will charged extra as Applicable. However of value above 10,000/- will be supplied on F.O.R. destination
Complaints	:	All complaints & enquiries at the customers request at their expenses.
Claims	:	All claims subject to JAIPUR Jurisdiction only.



PRECISION BIOMED PRIVATE LIMITED

PLOT NO. 193, SILVER SOIL INDUSTRIAL PARK,
VILLAGE : ANATPURA – CHIMANPURA, TEHSIL – CHOMU,
JAIPUR – 303702, RAJASTHAN (INDIA)

PH. - +91-9982806050 **CUSTOMER CARE : 07820806050**

Email : info@precisionbiomed.in Web : www.precisionbiomed.in

